

Apollo Anesthesia Workstation Daily Pre-Use Check-out Protocol

(Using On-Screen Self-Test)

1. Prepare machine for pre-use check
 - a. If the machine is on, turn it off
 - b. If gas was flowing overnight, change the CO₂ absorbent
2. Turn machine on
3. Perform manual check-out steps listed on screen
4. Open vaporizer of preferred agent to lowest marked concentration value and start automated check for:
 - a. Leak/compliance tests
 - b. Electronics tests
 - c. Agent analyzer tests
 - d. Flow sensor and piston checks
5. Test vaporizers
 - a. Test exclusion system – ensure that only one vaporizer can be turned on at a time
 - b. If desired, run leak test from ‘STANDBY’ screen with other vaporizer(s) set to lowest concentration
6. Verify operation of unidirectional valves, APL, and O₂ Flush
 - a. * Occlude y-piece
 - b. * Close APL and turn off all gas flows
 - c. * Use O₂ flush to fill breathing bag and generate positive circuit pressure
 - d. Turn on O₂ + other gas flows, check for smooth operation of flowmeter and verify increasing circuit pressure
- e. Slowly open APL valve and verify that pressure decreases smoothly
- f. Test interlink between O₂ & N₂O controls that protect against hypoxic O₂/N₂O mixtures
7. * Prepare machine for use
 - a. Set all vaporizers to zero
 - b. Turn off all flows and verify APL valve is open
 - c. Set up breathing circuit for use
8. * Check for necessary monitoring
 - a. Ensure all monitors are in place and functional

Abbreviated Pre-Use Check for Performance Between Anesthesia Cases (assumes machine has not been disturbed, for example by moving to another operating room)

1. Perform steps with a *
2. Replace patient suction, if necessary
3. Perform leak and compliance test from Standby screen, if necessary

Note: Developed for Apollo software version 3.21
Massachusetts General Hospital, Boston, MA
October 9, 2007